USERS' PROGRAM LIBRARY EUROPE
EUROPÄISCHE BENÜTZER-PROGRAMMBIBLIOTHEK
BIBLIOTHÈQUE EUROPÉENNE DE PROGRAMMES UTILISATEURS
LIBRERIA EUROPEA DEGLI UTILIZZATORI

Page 1 of 2

HP 41 PROGRAM SUBMITTAL FORM
PROGRAMMFORMBLATT/DOCUMENTATION DU PROGRAMME/GENERALITÀ SUL PROGRAMMA

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Operating limits and Warnings

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Limites et restrictions

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If you want to change a diameter or a % after pressing E, store 1 ST+16.

n is lowered by one in line 045 and 046 and will

be one to low when this is not corrected.

This program has been verified only with respect to the numerical example given in Program Description II. User accepts and uses this program material AT HIS OWN RISK, in reliance solely upon his own inspection of the program material and without reliance upon any representation or description concerning the program material.

NEITHER HP NOR THE CONTRIBUTOR MAKES ANY EXPRESS OR IMPLIED WARRANTY OF ANY KIND WITH REGARD TO THIS PROGRAM MATERIAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NEITHER HP NOR THE CONTRIBUTOR SHALL BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING OUT OF THE FURNISHING, USE OR PERFORMANCE OF THIS PROGRAM MATERIAL.

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Page 3 of 8

Example / Beispiel / Esempio: The R.I.D. plate is made in 2% agar noble in 0,9% NaCl.
The agar is melted in an boiling waterbath. After dissolving it is cooled till 56°C, the antiserum is added and after mixing it is poured between two glassplates from which one is siliconated.
When the agar is congealed, holes of approx. 2,5 mm are made in the agar and filled with antigen dilution. 100% is 1+0, 50% is 1+1 a.s.o.
After 24-72 hours (the 100% ring should be between 5 and 7 mm) the ring diameter is measured.
Because diameter² vs concentration has a linear relation a straight line can be calculated by linear regression, with $t_{95} = a$ and a correlation R. Then the titer is calculated.

For the example we made a R.I.D. plate with 0,150 ml antiserum in 12 ml agar, with a thickness of 0,15 cm.

We filled the holes with 3 µl antigen (1 mg/ml) and we found the next day the following ring diameters:

sketch
skizze
Croquis
Schizzo

dilution		diameter in mm
1+0	100%	7,0
1+1	50%	5,3
1+3	25%	4,0
1+9	10%	3,0

(Example continued)

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(Exemple, suite)

(Esempio)

PROGRAM DESCRIPTION III

 PROGRAMMBESCHREIBUNG III
 DESCRIPTION DU PROGRAMME III
 DESCRIZIONE DEL PROGRAMMA III

Page 7 of 8

Data Input Dateneingabe Données Dati	Keystrokes Tastendruck Touches Tasti	Result Résultat Résultat Risultato	Comments Kommentare Commentaires Commenti
	XEQ "TI"	UL AG?	3 µl (3.10 ⁻⁶ l) antigen was used
3	R/S	(AG) MG/ML	The antigen concentration of the 100% was 1 mg/ml
1	R/S	DIAM Y1	The diameter of the 100% ring is 7,0 mm
7	R/S	% X1	Display as input control
100	R/S	" 7,0/ 100%"	
		DIAM Y2	
5,3	R/S	% X2	
50	R/S	DIAM Y3	
4	R/S	% X3	
25	R/S	DIAM Y4	
3	R/S	% X4	
10	R/S	DIAM Y5	
			Press E (user mode !) to start the calculation of R, a and b
	E (LN)	R=0,9994	Correlation
	R/S	a=147,99	Tge
		b=4,99	Intercept.SQRT b should be the diameter of an empty hole (this is also a control!)
	R/S	ML AB?	0,150 ml antibody was used
0,15	R/S	ML AGAR?	This was mixed with 12 ml agar
12	R/S	MM AGAR	The agar layer was 0,15 cm thick
1,5	R/S		
		TITER=4,65	The titer is 4,65 mg/ml
			So one ml antiserum can bind 4,65 mg antigen.
			Press C to restart the program at the point R=0,9994
	A (ε+)	DIAM Y?	Press A when something goes wrong during data input
False diameter	R/S	%?	
False %	R/S	" Y/ X"	
		DIAM YA	Correction can also be done after pressing E.
	R/S	% X4	Do not forget to store 1 ST+16 before pressing E again.

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