

BORDER-LOCK® IMPRESSION TRAYS

THE FIRST SEMI CUSTOM TRAY

Developed by Schreinemakers

Disposable Border-Lock® impression trays: Dentulous & Edentulous



Normal

Square

Gothic



- Development of optimal Dynamic Pressure
- Uses 40% less impression material
- Loss of expensive metal trays is eliminated

- Fit identical to original Schreinemakers metal trays
- Excellent stability- rivaling that of metal trays
- Semi custom shape: Normal, Square, Gothic

Autoclavable Border-Lock® impression trays: Dentulous & Edentulous



High palate and postdam



Easy tray selection



Suction type retention guaranteed

- Fit identical to original Schreinemakers metal trays
- Tray can be heated up to 260° Celcius
- Excellent stability- rivaling that of metal trays

- Uses 40% less impression material – reducing costs
- Corrosion free material with 2 year Guarantee
- Excellent mechanical retention

Implant Border-Lock® impression trays: Dentulous & Edentulous



Position of abutment



Drill hole in impression tray



Final impression

- Fit identical to original Schreinemakers metal trays
- The drill holes can be made immediately at chairside
- The position of the abutment is perfectly visible

- Always a fitting tray available
- Tray is transparent
- Excellent control of costs and stock

Options in respect of Border-Lock® impression tray orders:

DISPOSABLE BORDER-LOCK IMPRESSION TRAYS **Introductory pack 96 border-lock trays dentulous** (1x set 8 upper + 8 lower + refills 10x nr 42-45, 51-54)

Art.No.

5401+5402



Storage Box 8 upper + Storage Box 8 lower dentulous

5411+5412

Refill box: 10 dentulous trays per size

5442-49, 5451-58



Introductory pack 40 border-lock trays edentulous (1x set upper + 5 lower + refills 3x nr 21-25, 31-35)

6400

Storage Box 5 upper + 5 lower edentulous

6410

Refill box: 5 edentulous trays per size

6421-25, 6431-35



AUTOCLAVABLE BORDER-LOCK IMPRESSION TRAYS **Set 8 upper + 8 lower border-lock trays dentulous**

5800

Single tray dentulous

5842-49, 5851-58



Set 5 upper + 5 lower border-lock trays edentulous Single tray edentulous

6800

6821-25, 6831-35



IMPLANT BORDER-LOCK IMPRESSION TRAYS

Introductory pack 96 implant trays dentulous (1x set 8 upper + 8 lower + refills 10x nr 42-45, 51-54)

5601+5602

Storage Box 8 upper + Storage Box 8 lower dentulous

5611+5612

Refill box: 10 dentulous trays per size

5642-49, 5651-58



Introductory pack 40 implant trays edentulous (1x set 5 upper + 5 lower + refills 3x nr 21-25, 31-35)

6600

Storage Box 5 upper + 5 lower edentulous

6610

Refill box: 5 edentulous trays per size

6621-25, 6631-35

